

Product Name: MST4 (9N11) Rabbit Monoclonal Antibody
Catalog #: AMRe14188

Summary

Production Name	MST4 (9N11) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	STK26
Alternative Names	MASK; Mst4; STE20 like kinase MST4; STK26;
Gene ID	51765.0
SwissProt ID	Q9P289.A synthetic peptide of human MST4

Application

Dilution Ratio	WB: 1:5000-1:20000
Molecular Weight	47kDa


Background

Product Name: MST4 (9N11) Rabbit Monoclonal Antibody
Catalog #: AMRe14188

Mediator of cell growth. Modulates apoptosis. Mediator of cell growth (PubMed:11641781, PubMed:17360971). Modulates apoptosis (PubMed:11641781, PubMed:17360971). In association with STK24 negatively regulates Golgi reorientation in polarized cell migration upon RHO activation (PubMed:27807006).

Research Area

Image Data



Western blot detection of MST4 in K562,C6 cell lysates using MST4 antibody(1:1000 diluted).

Note

For research use only.