

Product Name: Nogo B Receptor (19G15) Rabbit Monoclonal Antibody
Catalog #: AMRe14785

Summary

Production Name	Nogo B Receptor (19G15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	NUS1
Alternative Names	
Gene ID	116150.0
SwissProt ID	Q96E22.A synthetic peptide of human Nogo B receptor

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	33kDa


Background

**Product Name: Nogo B Receptor (19G15) Rabbit
Monoclonal Antibody
Catalog #: AMRe14785**



Research Area

Image Data



Western blot detection of Nogo B receptor in K562,3T3,Hela cell lysates using Nogo B receptor antibody(1:1000 diluted).

Note

For research use only.