

Product Name: Oct-2 (14K1) Rabbit Monoclonal Antibody
Catalog #: AMRe15100



Summary

Production Name	Oct-2 (14K1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	POU2F2
Alternative Names	
Gene ID	5452.0
SwissProt ID	P09086.A synthetic peptide of human Oct-2

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	51kDa


Background

Product Name: Oct-2 (14K1) Rabbit Monoclonal Antibody
Catalog #: AMRe15100



Research Area

Image Data



Western blot detection of Oct-2 in A549,HL-60,U2OS cell lysates using Oct-2 antibody(1:1000 diluted).

Note

For research use only.