

Product Name: p53 (9D15) Rabbit Monoclonal Antibody
Catalog #: AMRe15641



Summary

Production Name	p53 (9D15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	Tp53
Alternative Names	
Gene ID	22059.0
SwissProt ID	P02340.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	43 kDa


Background

Product Name: p53 (9D15) Rabbit Monoclonal Antibody
Catalog #: AMRe15641



Research Area

Image Data



Western blot analysis of p53 in 3T3 lysates using p53 antibody.

Note

For research use only.