

Product Name: PAX8 (2K7) Rabbit Monoclonal Antibody
Catalog #: AMRe15795



Summary

Production Name	PAX8 (2K7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PAX8
Alternative Names	Paired box 8; Paired box gene 8; paired box homeotic gene 8; Paired box protein Pax 8; Paired domain gene 8; PAX8;
Gene ID	7849.0
SwissProt ID	Q06710.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	48kDa

Product Name: PAX8 (2K7) Rabbit Monoclonal Antibody
Catalog #: AMRe15795




Background

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells. Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

Research Area

Image Data



Western blot analysis of PAX8 expression in SK-OV-3 cell lysate.

Note

For research use only.