

Product Name: POLG (9R12) Rabbit Monoclonal Antibody
Catalog #: AMRe16350

Summary

Production Name	POLG (9R12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	POLG
Alternative Names	
Gene ID	5428.0
SwissProt ID	P54098.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	140 kDa


Background

Product Name: POLG (9R12) Rabbit Monoclonal Antibody
Catalog #: AMRe16350



Research Area

Image Data



Western blot analysis of POLG in HeLa, A549, HL-60, U2OS lysates using POLG antibody.

Note

For research use only.