

Product Name: PPP2R5E (10K3) Rabbit Monoclonal Antibody
Catalog #: AMRe16431

Summary

Production Name	PPP2R5E (10K3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PPP2R5E
Alternative Names	PPP2R5E;
Gene ID	5529.0
SwissProt ID	Q16537.A synthetic peptide of human PPP2R5E

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	55kDa

Background


Product Name: PPP2R5E (10K3) Rabbit Monoclonal Antibody
Catalog #: AMRe16431



The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment. The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

Research Area

Image Data



Western blot detection of PPP2R5E in A549,HL-60,U251 cell lysates using PPP2R5E antibody(1:1000 diluted).

Note

For research use only.