

Product Name: PWP1 (16C6) Rabbit Monoclonal Antibody
Catalog #: AMRe16703

Summary

Production Name	PWP1 (16C6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PWP1
Alternative Names	PWP1;
Gene ID	11137.0
SwissProt ID	Q13610.A synthetic peptide of human PWP1

Application

Dilution Ratio	WB: 1:1000-1:2000
Molecular Weight	56kDa


Background

Product Name: PWP1 (16C6) Rabbit Monoclonal Antibody
Catalog #: AMRe16703

May play an important role in cell growth and/or transcription. Chromatin-associated factor that regulates transcription (PubMed:29065309). Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription (PubMed:29065309). Regulates the epigenetic status of rDNA (PubMed:29065309).

Research Area

Image Data



Western blot detection of PWP1 in K562 cell lysates using PWP1 antibody(1:1000 diluted).

Note

For research use only.