

Product Name: RAD52 (13D10) Rabbit Monoclonal Antibody
Catalog #: AMRe16847

Summary

Production Name	RAD52 (13D10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	RAD52
Alternative Names	Rad52; RAD52 homolog;
Gene ID	5893.0
SwissProt ID	P43351.A synthetic peptide of human RAD52

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	46kDa


Background

Product Name: RAD52 (13D10) Rabbit Monoclonal Antibody
Catalog #: AMRe16847

Involved in double-stranded break repair. Plays a central role in genetic recombination and DNA repair by promoting the annealing of complementary single-stranded DNA and by stimulation of the RAD51 recombinase. Involved in double-stranded break repair. Plays a central role in genetic recombination and DNA repair by promoting the annealing of complementary single-stranded DNA and by stimulation of the RAD51 recombinase.

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using RM6417 at 1:1000.

Note

For research use only.