

Product Name: RhoGAP (11C13) Rabbit Monoclonal Antibody
Catalog #: AMRe17131



Summary

Production Name	RhoGAP (11C13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ARHGAP5
Alternative Names	ARHGAP5; p190 B; p190 RhoGAP; Rho GTPase activating protein 5; RhoGAP; RhoGAP5;
Gene ID	394.0
SwissProt ID	Q13017.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	172kDa

Product Name: RhoGAP (11C13) Rabbit Monoclonal Antibody
Catalog #: AMRe17131




Background

GTPase-activating protein for Rho family members (PubMed:8537347). GTPase-activating protein for Rho family members (PubMed:8537347).

Research Area

Image Data



Western blot analysis of RhoGAP expression in HeLa cell lysate.

Note

For research use only.