

Product Name: S100 alpha9 (13Z5) Rabbit Monoclonal Antibody
Catalog #: AMRe17462



Summary

Production Name	S100 alpha9 (13Z5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	S100A9
Alternative Names	
Gene ID	20202.0
SwissProt ID	P31725.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	13 kDa


Background

Product Name: S100 alpha9 (13Z5) Rabbit Monoclonal Antibody
Catalog #: AMRe17462



Research Area

Image Data



Western blot analysis of S100A9 in mouse spleen, mouse lung lysates using S100 alpha9 antibody.

Note

For research use only.