

**Product Name: Salivary alpha amylase (5D14) Rabbit
Monoclonal Antibody
Catalog #: AMRe17579**



Summary

Production Name	Salivary alpha amylase (5D14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	Alpha amylase; AMY1; AMY1C; AMY2; AMY2B; AMY3; HXA;
Gene ID	
SwissProt ID	P04745.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	58kDa

**Product Name: Salivary alpha amylase (5D14) Rabbit
Monoclonal Antibody
Catalog #: AMRe17579**




Background

Belongs to the glycosyl hydrolase 13 family.

Research Area

Image Data



Western blot analysis of Salivary alpha amylase expression in Human pancreas lysate.

Note

For research use only.