

Product Name: SEC61A (18P13) Rabbit Monoclonal Antibody
Catalog #: AMRe17694

Summary

Production Name	SEC61A (18P13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	fb62c11; HSEC61; SEC 61; SEC 61A; Sec61 alpha 1; SEC61; SEC61A1; SEC61A2;
Gene ID	
SwissProt ID	P61619/Q9H9S3.Recombinant protein of human SEC61A

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	52kDa

Background


Product Name: SEC61A (18P13) Rabbit Monoclonal Antibody
Catalog #: AMRe17694



Plays a crucial role in the insertion of secretory and membrane polypeptides into the ER. Required for assembly of membrane and secretory proteins. Tightly associated with membrane-bound ribosomes, either directly or through adapter proteins.

Research Area

Image Data



Western blot detection of SEC61A in K562,C6,3T3 cell lysates using SEC61A antibody(1:1000 diluted).

Note

For research use only.