

**Product Name: SH2D1A (7N13) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe17834**

---

## Summary

|                        |                                          |
|------------------------|------------------------------------------|
| <b>Production Name</b> | SH2D1A (7N13) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal Antibody               |
| <b>Host</b>            | Rabbit                                   |
| <b>Application</b>     | WB                                       |
| <b>Reactivity</b>      | Human,Mouse,Rat                          |

## Performance

|                     |                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                                 |
| <b>Modification</b> | Unmodified                                                                                                   |
| <b>Isotype</b>      | IgG                                                                                                          |
| <b>Clonality</b>    | Monoclonal                                                                                                   |
| <b>Form</b>         | Liquid                                                                                                       |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                     |
| <b>Buffer</b>       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA. |
| <b>Purification</b> | Affinity purification                                                                                        |

## Immunogen

|                          |                                                |
|--------------------------|------------------------------------------------|
| <b>Gene Name</b>         | SH2D1A                                         |
| <b>Alternative Names</b> |                                                |
| <b>Gene ID</b>           | 4068.0                                         |
| <b>SwissProt ID</b>      | O60880.A synthetic peptide of human SH2D1A/SAP |

## Application

|                         |                   |
|-------------------------|-------------------|
| <b>Dilution Ratio</b>   | WB: 1:1000-1:2000 |
| <b>Molecular Weight</b> | 14kDa             |

## Background

---


**Product Name: SH2D1A (7N13) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe17834**

---



## Research Area

## Image Data



Western blot detection of SH2D1A/SAP in MCF7,C6 cell lysates using SH2D1A/SAP antibody(1:1000 diluted).

## Note

For research use only.