

Product Name: STIP1 (16V12) Rabbit Monoclonal Antibody
Catalog #: AMRe18387

Summary

Production Name	STIP1 (16V12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	STIP1
Alternative Names	Hop; P60; STI1; STI1L; STIP1; Stress induced phosphoprotein 1;
Gene ID	10963.0
SwissProt ID	P31948.A synthetic peptide of human STIP1

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	63kDa


Background

Product Name: STIP1 (16V12) Rabbit Monoclonal Antibody
Catalog #: AMRe18387

Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB). Acts as a co-chaperone for HSP90AA1 (PubMed: <http://www.uniprot.org/citations/27353360>). Mediates the association of the molecular chaperones HSPA8/HSC70 and HSP90 (By similarity).

Research Area

Image Data



Western blot detection of STIP1 in K562,Rat Brain,C6 cell lysates using STIP1 antibody(1:1000 diluted).

Note

For research use only.