

**Product Name: Superoxide Dismutase 1 (10Y2) Rabbit
Monoclonal Antibody
Catalog #: AMRe18447**



Summary

Production Name	Superoxide Dismutase 1 (10Y2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	SOD1
Alternative Names	
Gene ID	6647.0
SwissProt ID	P00441.A synthetic peptide of human Superoxide Dismutase 1

Application

Dilution Ratio	WB: 1:5000
Molecular Weight	16kDa


Background

**Product Name: Superoxide Dismutase 1 (10Y2) Rabbit
Monoclonal Antibody
Catalog #: AMRe18447**



Research Area

Image Data



Western blot detection of Superoxide Dismutase 1 in K562,Rat Brain,C6,3T3,Hela cell lysates using Superoxide Dismutase 1 antibody(1:1000 diluted).

Note

For research use only.