

Product Name: Syntaxin 16 (10E13) Rabbit Monoclonal Antibody
Catalog #: AMRe18512



Summary

Production Name	Syntaxin 16 (10E13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	STX16
Alternative Names	hsyn16; Stx16; Syn16;
Gene ID	8675.0
SwissProt ID	O14662.A synthetic peptide of human Syntaxin 16

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	37kDa

Background


Product Name: Syntaxin 16 (10E13) Rabbit Monoclonal Antibody
Catalog #: AMRe18512



SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network. SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.

Research Area

Image Data



Western blot detection of Syntaxin 16 in K562,Rat Brain lysates using Syntaxin 16 antibody(1:1000 diluted).

Note

For research use only.