

Product Name: Syntaxin 3 (18G14) Rabbit Monoclonal Antibody
Catalog #: AMRe18513



Summary

Production Name	Syntaxin 3 (18G14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	STX3
Alternative Names	Stx3; STX3A; Syntaxin 3A;
Gene ID	6809.0
SwissProt ID	Q13277.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	33kDa

Product Name: Syntaxin 3 (18G14) Rabbit Monoclonal Antibody
Catalog #: AMRe18513




Background

Potentially involved in docking of synaptic vesicles at presynaptic active zones. Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Research Area

Image Data



Western blot analysis of Syntaxin 3 expression in HepG2 cell lysate.

Note

For research use only.