

Product Name: TMEM192 (10L14) Rabbit Monoclonal Antibody
Catalog #: AMRe19054



Summary

Production Name	TMEM192 (10L14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TMEM192
Alternative Names	TMEM192; Transmembrane protein 192;
Gene ID	201931.0
SwissProt ID	Q8IY95.Recombinant protein of human TMEM192

Application

Dilution Ratio	WB: 1:1000-1:2000
Molecular Weight	31kDa

Background


Product Name: TMEM192 (10L14) Rabbit Monoclonal Antibody
Catalog #: AMRe19054



Belongs to the TMEM192 family.

Research Area

Image Data



Western blot detection of TMEM192 in K562 cell lysates using TMEM192 antibody(1:1000 diluted).

Note

For research use only.