

**Product Name: Transcriptional Enhancer Factor TEF
(14G4) Rabbit Monoclonal Antibody**
Catalog #: AMRe19199

Summary

Production Name	Transcriptional Enhancer Factor TEF (14G4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TEAD4
Alternative Names	
Gene ID	7004.0
SwissProt ID	Q15561.Recombinant protein of human TEAD4

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	49kDa


Background

**Product Name: Transcriptional Enhancer Factor TEF
(14G4) Rabbit Monoclonal Antibody**
Catalog #: AMRe19199



Research Area

Image Data



Western blot detection of TEAD1/TEAD2/TEAD3/TEAD4 in K562, HeLa cell lysates using TEAD1/TEAD2/TEAD3/TEAD4 antibody (1:1000 diluted).

Note

For research use only.