

Product Name: Transmembrane Protein 161B (19D14)
Rabbit Monoclonal Antibody
Catalog #: AMRe19219



Summary

Production Name	Transmembrane Protein 161B (19D14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	TMEM161B
Alternative Names	
Gene ID	153396.0
SwissProt ID	Q8NDZ6.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	55 kDa


Background

Product Name: Transmembrane Protein 161B (19D14)
Rabbit Monoclonal Antibody
Catalog #: AMRe19219



Research Area

Image Data



Western blot analysis of TMEM161B in K562, HeLa lysates using Transmembrane Protein 161B antibody.

Note

For research use only.