

Product Name: Transmembrane Protein 43 (9N13)
Rabbit Monoclonal Antibody
Catalog #: AMRe19220



Summary

Production Name	Transmembrane Protein 43 (9N13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TMEM43
Alternative Names	
Gene ID	79188.0
SwissProt ID	Q9BTV4.A synthetic peptide of human TMEM43

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	45kDa


Background

Product Name: Transmembrane Protein 43 (9N13)
Rabbit Monoclonal Antibody
Catalog #: AMRe19220



Research Area

Image Data



Western blot detection of TMEM43 in K562,C6 cell lysates using TMEM43 antibody(1:1000 diluted).

Note

For research use only.