

Product Name: TWEAKR (14X18) Rabbit Monoclonal Antibody
Catalog #: AMRe19442



Summary

Production Name	TWEAKR (14X18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	TNFRSF12A
Alternative Names	CD266 antigen; FGF inducible 14; Fibroblast growth factor inducible immediate early response protein 14; FN14; TNFRSF 12A; Tweak receptor; TweakR;
Gene ID	51330.0
SwissProt ID	Q9NP84.

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	14kDa

Product Name: TWEAKR (14X18) Rabbit Monoclonal Antibody
Catalog #: AMRe19442




Background

Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins. Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins.

Research Area

Image Data



Western blot analysis of TWEAKR expression in HUVEC cell lysate.

Note

For research use only.