

Product Name: VEGFD (7C13) Rabbit Monoclonal Antibody
Catalog #: AMRe19774

Summary

Production Name	VEGFD (7C13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	VEGFD
Alternative Names	FIGF; VEGFD;
Gene ID	2277.0
SwissProt ID	O43915.A synthetic peptide of human VEGFD

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	40kDa


Background

Product Name: VEGFD (7C13) Rabbit Monoclonal Antibody
Catalog #: AMRe19774

Growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. Growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. May function in the formation of the venous and lymphatic vascular systems during embryogenesis, and also in the maintenance of differentiated lymphatic endothelium in adults. Binds and activates VEGFR-2 (KDR/FLK1) and VEGFR-3 (FLT4) receptors.

Research Area

Image Data



Western blot detection of VEGFD in Rat Brain,C6,HeLa cell lysates using VEGFD antibody(1:1000 diluted).

Note

For research use only.