

Product Name: Versican (15C15) Rabbit Monoclonal Antibody
Catalog #: AMRe19779



Summary

Production Name	Versican (15C15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	VCAN
Alternative Names	CSPG2; ERVR; GHAP; PGM; VCAN; Versican; WGN1;
Gene ID	1462.0
SwissProt ID	P13611.Recombinant protein of human Versican

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	373kDa

Background


Product Name: Versican (15C15) Rabbit Monoclonal Antibody
Catalog #: AMRe19779



May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid. May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid.

Research Area

Image Data



Western blot analysis of extracts from Mouse brain tissue using RM5694 at 1:1000.

Note

For research use only.