

Product Name: Cyclooxygenase 1 Rabbit Polyclonal Antibody
Catalog #: APRab00091

Summary

Production Name	Cyclooxygenase 1 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	PTGS1
Alternative Names	PTGS1; COX1; Prostaglandin G/H synthase 1; Cyclooxygenase-1; COX-1; Prostaglandin H2 synthase 1; PGH synthase 1; PGHS-1; PHS 1; Prostaglandin-endoperoxide synthase 1
Gene ID	5742
SwissProt ID	P23219

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-
-----------------------	--

Product Name: Cyclooxygenase 1 Rabbit Polyclonal Antibody
Catalog #: APRab00091



1/100

Molecular Weight

Calculated MW: 69 kDa; Observed MW: 69 kDa


Background

Cyclooxygenase-1 May play an important role in regulating or promoting cell proliferation in some normal and neoplastically transformed cells. Belongs to the prostaglandin G/H synthase family. Homodimer. 2 isoforms of the human protein are produced by alternative splicing.


Research Area

Immunology

Image Data



Western blot analysis of COX1 in A431 lysates using Cyclooxygenase 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using COX1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

**Product Name: Cyclooxygenase 1 Rabbit Polyclonal
Antibody**
Catalog #: APRab00091



For research use only.