

Product Name: INPP4B Rabbit Polyclonal Antibody
Catalog #: APRab00101



Summary

Production Name	INPP4B Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	INPP4B 4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa;
Alternative Names	Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3
Gene ID	8821
SwissProt ID	O15327

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100
Molecular Weight	Calculated MW: 105 kDa; Observed MW: 105 kDa

Product Name: INPP4B Rabbit Polyclonal Antibody
Catalog #: APRab00101




Background

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

Research Area

Signal Transduction

Image Data



Western blot analysis of INPP4B in MCF-7 lysates using INPP4B antibody.

Note

For research use only.