

**Product Name: IL-1 Receptor Antagonist Protein Rabbit Polyclonal Antibody**  
**Catalog #: APRab00104**

---

## Summary

<b>Production Name</b>	IL-1 Receptor Antagonist Protein Rabbit Polyclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P,ICC/IF,FC,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification</b>	Affinity Chromatography

## Immunogen

<b>Gene Name</b>	IL1RN
<b>Alternative Names</b>	IL-1 Receptor Antagonist Protein; DIRA; IRAP; IL1F3; IL1RA; MVCD4; IL-1RN; IL-1ra; IL-1ra3; ICIL-1RA
<b>Gene ID</b>	3557
<b>SwissProt ID</b>	P18510

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
-----------------------	---

**Product Name: IL-1 Receptor Antagonist Protein Rabbit  
Polyclonal Antibody  
Catalog #: APRab00104**



---

**Molecular Weight**

Calculated MW: 20 kDa; Observed MW: 20 kDa


**Background**

Interleukin-1 (IL-1) is a key pro-inflammatory cytokine that has diverse actions in the brain as a regulator of host defense responses and a mediator of inflammation. Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity.

**Research Area**

Immunology

**Image Data**



Western blot analysis of IL-1RA in (1) HeLa lysates; (2) 3T3 lysates using IL-1 Receptor Antagonist Protein antibody.

**Note**

For research use only.