

Product Name: TCTP Rabbit Polyclonal Antibody
Catalog #: APRab00222



Summary

Production Name	TCTP Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	TPT1
Alternative Names	Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23
Gene ID	7178
SwissProt ID	P13693

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Product Name: TCTP Rabbit Polyclonal Antibody
Catalog #: APRab00222




Background

Involved in calcium binding and microtubule stabilization.

Research Area

Immunology

Image Data



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.

Note

For research use only.