

Product Name: Cyclophilin F Rabbit Polyclonal Antibody
Catalog #: APRab00232



Summary

Production Name	Cyclophilin F Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	PPIF
Alternative Names	PPIF; CYP3; Peptidyl-prolyl cis-trans isomerase F; mitochondrial; PPlase F; Cyclophilin F; Rotamase F
Gene ID	10105
SwissProt ID	P30405

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 18 kDa

Product Name: Cyclophilin F Rabbit Polyclonal Antibody
Catalog #: APRab00232




Background

PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Research Area

Immunology

Image Data



Western blot analysis of Cyclophilin F in HeLa lysates using Cyclophilin F antibody.

Note

For research use only.