

Product Name: Thymidylate Synthase Rabbit Polyclonal Antibody
Catalog #: APRab00246

Summary

Production Name	Thymidylate Synthase Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	TYMS
Alternative Names	TS; TMS; HST422
Gene ID	7298
SwissProt ID	P04818

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 36 kDa

Product Name: Thymidylate Synthase Rabbit Polyclonal Antibody
Catalog #: APRab00246


Background

Contributes to the de novo mitochondrial thymidylate biosynthesis pathway.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Thymidylate Synthase in HeLa lysates using Thymidylate Synthase antibody.

Note

For research use only.