

Product Name: Aspartate Aminotransferase Rabbit Polyclonal Antibody
Catalog #: APRab00263

Summary

Production Name	Aspartate Aminotransferase Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	GOT1
Alternative Names	Aspartate aminotransferase; cytoplasmic; Glutamate oxaloacetate transaminase 1; Transaminase A
Gene ID	2805
SwissProt ID	P17174

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa

**Product Name: Aspartate Aminotransferase Rabbit
Polyclonal Antibody**
Catalog #: APRab00263




Background

Plays a key role in amino acid metabolism.

Research Area

Signal Transduction

Image Data



Western blot analysis of Aspartate Aminotransferase in MCF-7 lysates using Aspartate Aminotransferase antibody.

Note

For research use only.