

**Product Name: Phospho-c-Jun (Ser63) Rabbit Polyclonal Antibody**  
**Catalog #: APRab00710**

---

## Summary

<b>Production Name</b>	Phospho-c-Jun (Ser63) Rabbit Polyclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P,ICC/IF
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phosphorylated
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification</b>	Affinity Chromatography

## Immunogen

<b>Gene Name</b>	JUN
<b>Alternative Names</b>	Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
<b>Gene ID</b>	3725
<b>SwissProt ID</b>	P05412

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
<b>Molecular Weight</b>	Calculated MW: 36 kDa; Observed MW: 48 kDa

**Product Name: Phospho-c-Jun (Ser63) Rabbit Polyclonal Antibody**  
**Catalog #: APRab00710**

---


## Background

c-jun is a transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'. Promotes activity of NR5A1 when phosphorylated by HIPK3 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation.


## Research Area

Epigenetics and Nuclear Signaling

## Image Data




Western blot analysis of c-Jun Phosphorylation in NIH/3T3 lysates treated with Anisomycin using Phospho-c-Jun (Ser63) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using Phospho-c-Jun (S63) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**Product Name: Phospho-c-Jun (Ser63) Rabbit Polyclonal Antibody**  
**Catalog #: APRab00710**

---



Immunofluorescence analysis of Phospho-c-Jun (Ser63) in HeLa treated with anisomycin using Phospho-c-Jun (S63) antibody.

### **Note**

For research use only.