

Product Name: Phospho-Rac1 (Ser71) Rabbit Polyclonal Antibody
Catalog #: APRab00917

Summary

Production Name	Phospho-Rac1 (Ser71) Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	RAC1
Alternative Names	Ras-related C3 botulinum toxin substrate 1; Cell migration-inducing gene 5 protein; Ras-like protein TC25; p21-Rac1; RAC1; TC25
Gene ID	5879
SwissProt ID	P63000

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 21 kDa; Observed MW: 28 kDa

Product Name: Phospho-Rac1 (Ser71) Rabbit Polyclonal Antibody
Catalog #: APRab00917


Background

Plasma membrane-associated small GTPase which cycles between active GTP-bound and inactive GDP-bound states.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phospho-Rac1 (Ser71) in 293 lysates using Phospho-Rac1 (Ser71) antibody.

Note

For research use only.