

Product Name: Thrombospondin 2 Rabbit Polyclonal Antibody
Catalog #: APRab03354



Summary

Production Name	Thrombospondin 2 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	THBS2
Alternative Names	THBS2; TSP2; Thrombospondin-2
Gene ID	7058
SwissProt ID	P35442

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000
Molecular Weight	Calculated MW: 130 kDa; Observed MW: 130 kDa

Background

Product Name: Thrombospondin 2 Rabbit Polyclonal Antibody
Catalog #: APRab03354




Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions.


Research Area

Cardiovascular

Image Data



Western blot analysis of Thrombospondin 2 in mouse heart lysates using Thrombospondin 2 antibody.



Immunohistochemistry analysis of paraffin-embedded rat testis using Thrombospondin 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.