

Product Name: PPAR alpha Rabbit Polyclonal Antibody
Catalog #: APRab03383



Summary

Production Name	PPAR alpha Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	PPARA
Alternative Names	PPARA; NR1C1; PPAR; Peroxisome proliferator-activated receptor alpha; PPAR-alpha; Nuclear receptor subfamily 1 group C member 1
Gene ID	5465
SwissProt ID	Q07869

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

Product Name: PPAR alpha Rabbit Polyclonal Antibody
Catalog #: APRab03383




Ligand-activated transcription factor.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of PPAR alpha in 3T3 lysates using PPAR alpha antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen tissue using PPARα antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of PPAR alpha in HeLa cells using PPARα antibody (red) and Galectin3 antibody (green)

Product Name: PPAR alpha Rabbit Polyclonal Antibody
Catalog #: APRab03383



Note

For research use only.