

Product Name: FYN Rabbit Polyclonal Antibody
Catalog #: APRab03633



Summary

Production Name	FYN Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	FYN
Alternative Names	FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn
Gene ID	2534
SwissProt ID	P06241

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 61 kDa; Observed MW: 61 kDa

Background

Product Name: FYN Rabbit Polyclonal Antibody
Catalog #: APRab03633




Tyrosine-protein kinase implicated in the control of cell growth. Plays a role in the regulation of intracellular calcium levels, with isoform 2 showing the greater ability to mobilize cytoplasmic calcium in comparison to isoform 1. Required in brain development and mature brain function with important roles in the regulation of axon growth, axon guidance, and neurite extension.

Research Area

Neuroscience

Image Data



Western blot analysis of FYN in Ramos, Jurkat, C6, 3T3 and HeLa lysates using Fyn antibody.

Note

For research use only.