

Product Name: IL-1 beta Rabbit Polyclonal Antibody
Catalog #: APRab03854



Summary

Production Name	IL-1 beta Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	IL1B
Alternative Names	IL1B; IL1F2; Interleukin-1 beta; IL-1 beta; Catabolin
Gene ID	3553
SwissProt ID	P01584

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 31 kDa; Observed MW: 31 kDa

Background

Produced by activated macrophages, IL-1 stimulates thymocyte proliferation by inducing IL-2 release, B-cell maturation and

Product Name: IL-1 beta Rabbit Polyclonal Antibody
Catalog #: APRab03854




proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release of prostaglandin and collagenase from synovial cells.


Research Area

Immunology

Image Data



Western blot analysis of IL-1 beta in various cell lines lysates using IL-1B antibody.



Western blot analysis of IL-1 beta in various cell lines lysates using IL-1B antibody.

Note

For research use only.