

Product Name: Jun B (phospho Ser79) Rabbit Polyclonal Antibody
Catalog #: APRab04912

Summary

Production Name	Jun B (phospho Ser79) Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	JUNB
Alternative Names	JUNB; Transcription factor jun-B
Gene ID	3726.0
SwissProt ID	P17275.The antiserum was produced against synthesized peptide derived from human JunB around the phosphorylation site of Ser79. AA range:46-95

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:20000
Molecular Weight	


Background

Product Name: Jun B (phospho Ser79) Rabbit Polyclonal Antibody
Catalog #: APRab04912

function:Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3',induction:By growth factors.,similarity:Belongs to the bZIP family.,similarity:Belongs to the bZIP family. Jun subfamily.,similarity:Contains 1 bZIP domain.,subunit:Binds DNA as an homodimer or as an heterodimer with another member of the jun/fos family.,function:Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3',induction:By growth factors.,similarity:Belongs to the bZIP family.,similarity:Belongs to the bZIP family. Jun subfamily.,similarity:Contains 1 bZIP domain.,subunit:Binds DNA as an homodimer or as an heterodimer with another member of the jun/fos family.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using JunB (Phospho-Ser79) Antibody.
The picture on the right is blocked with the phospho peptide.

Note

For research use only.