

Product Name: STAT5A/B (Phospho-Ser725/730) Rabbit Polyclonal Antibody
Catalog #: APRab06088

Summary

Production Name	STAT5A/B (Phospho-Ser725/730) Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	STAT5A STAT5
Alternative Names	Signal transducer and activator of transcription 5A
Gene ID	6776.0
SwissProt ID	P42229.

Application


Dilution Ratio	WB 1:500-2000
Molecular Weight	87kD

Background

**Product Name: STAT5A/B (Phospho-Ser725/730) Rabbit
Polyclonal Antibody**
Catalog #: APRab06088

Research Area

Image Data



Western Blot analysis of heLa 293T using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10000 25°C, 1.5 hours

Note

For research use only.