

Product Name: ALKB7 Rabbit Polyclonal Antibody
Catalog #: APRab06785



Summary

Production Name	ALKB7 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ALKBH7 ABH7 SPATA11 UNQ6002/PRO34564
Alternative Names	
Gene ID	84266.0
SwissProt ID	Q9BT30.Synthesized peptide derived from human protein . at AA range: 10-90

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	24kD

Background


similarity:Belongs to the alkB family.,similarity:Belongs to the alkB family.,

Product Name: ALKB7 Rabbit Polyclonal Antibody
Catalog #: APRab06785



Research Area

Image Data



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

Note

For research use only.