

Product Name: ALKBH5 Rabbit Polyclonal Antibody
Catalog #: APRab06788



Summary

Production Name	ALKBH5 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ALKBH5 ABH5 OFOXD1
Alternative Names	Probable alpha-ketoglutarate-dependent dioxygenase ABH5 (EC 1.14.11.-;Alkylated DNA repair protein alkB homolog 5)
Gene ID	54890.0
SwissProt ID	Q6P6C2.Synthesized peptide derived from human ALKBH5

Application

Dilution Ratio	IHC 1:50-200
Molecular Weight	51kD

Background


similarity:Belongs to the alkB family.,similarity:Belongs to the alkB family.,

Product Name: ALKBH5 Rabbit Polyclonal Antibody
Catalog #: APRab06788



Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

Note

For research use only.