

**Product Name: Annexin VII Rabbit Polyclonal Antibody**  
**Catalog #: APRab06928**



## Summary

<b>Production Name</b>	Annexin VII Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ELISA
<b>Reactivity</b>	Human,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	ANXA7
<b>Alternative Names</b>	ANXA7; ANX7; SNX; OK/SW-cl.95; Annexin A7; Annexin VII; Annexin-7; Synexin
<b>Gene ID</b>	310.0
<b>SwissProt ID</b>	P20073.The antiserum was produced against synthesized peptide derived from the Internal region of human ANXA7. AA range:271-320

## Application

<b>Dilution Ratio</b>	WB 1:500 - 1:2000. IHC: 1:100-1:300. ELISA: 1:20000..
<b>Molecular Weight</b>	53kD

## Background


**Product Name: Annexin VII Rabbit Polyclonal Antibody**  
**Catalog #: APRab06928**




Annexin VII is a member of the annexin family of calcium-dependent phospholipid binding proteins. The Annexin VII gene contains 14 exons and spans approximately 34 kb of DNA. An alternatively spliced cassette exon results in two mRNA transcripts of 2.0 and 2.4 kb which are predicted to generate two protein isoforms differing in their N-terminal domain. The alternative splicing event is tissue specific and the mRNA containing the cassette exon is prevalent in brain, heart and skeletal muscle. The transcripts also differ in their 3'-non coding regions by the use of two alternative poly(A) signals. Annexin VII encodes a protein with a molecular weight of approximately 51 kDa with a unique, highly hydrophobic N-terminal domain of 167 amino acids and a conserved C-terminal region of 299 amino acids. The latter domain is composed of alternating hydrophobic and hydrophilic segments. Structural analysis of the domain: A pair of annexin repeats may form one binding site for calcium and phospholipid. function: Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis. similarity: Belongs to the annexin family. similarity: Contains 4 annexin repeats. tissue specificity: Isoform 1 is expressed in brain, heart and skeletal muscle. Isoform 2 is more abundant in liver, lung, kidney, spleen, fibroblasts and placenta.

## Research Area

## Image Data




Western blot analysis of lysate from A549 cells, using ANXA7 Antibody.




Western Blot analysis of A549 cells using Annexin VII Polyclonal Antibody.. Secondary antibody was diluted at 1:20000


**Product Name: Annexin VII Rabbit Polyclonal Antibody**  
**Catalog #: APRab06928**



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:100

## **Note**

For research use only.