

Product Name: Annexin X Rabbit Polyclonal Antibody
Catalog #: APRab06929



Summary

Production Name	Annexin X Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ANXA10
Alternative Names	ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14
Gene ID	11199.0
SwissProt ID	Q9UJ72.The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:20000
Molecular Weight	38kD

Background


Product Name: Annexin X Rabbit Polyclonal Antibody
Catalog #: APRab06929




This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],similarity:Belongs to the annexin family.,similarity:Contains 4 annexin repeats,

Research Area

Image Data



Western blot analysis of the lysates from Jurkat cells using ANXA10 antibody.



Western Blot analysis of various cells using Annexin X Polyclonal Antibody

Note

For research use only.