

Product Name: ARP19 Rabbit Polyclonal Antibody
Catalog #: APRab07157



Summary

Production Name	ARP19 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ARPP19
Alternative Names	
Gene ID	10776.0
SwissProt ID	P56211.Synthesized peptide derived from part region of human protein AA range: 1-50

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	12kD

Background

The 19-kD cAMP-regulated phosphoprotein plays a role in regulating mitosis by inhibiting protein phosphatase-2A (PP2A; see MIM 176915) (summary by Gharbi-Ayachi et al., 2010 [PubMed 21164014]).[supplied by OMIM, Feb 2011],function:May


Product Name: ARP19 Rabbit Polyclonal Antibody
Catalog #: APRab07157



enhance GAP-43 expression by binding to the NGF-regulatory region of its mRNA.,PTM:Phosphorylated by PKA.,similarity:Belongs to the endosulfine family.,

Research Area

Image Data



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

Note

For research use only.