

Product Name: BTF Rabbit Polyclonal Antibody
Catalog #: APRab07685



Summary

Production Name	BTF Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	BCLAF1
Alternative Names	BCLAF1; BTF; KIAA0164; Bcl-2-associated transcription factor 1; Btf
Gene ID	9774.0
SwissProt ID	Q9NYF8.The antiserum was produced against synthesized peptide derived from human BCLAF1. AA range:281-330

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:10000. Not yet tested in other applications.
Molecular Weight	106kD

Background


Product Name: BTF Rabbit Polyclonal Antibody
Catalog #: APRab07685




This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],function:Death-promoting transcriptional repressor.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with Bcl-2 related proteins, EMD, with the adenovirus E1B 19 kDa protein and with DNA.,tissue specificity:Ubiquitous.,

Research Area

Image Data




Western blot analysis of the lysates from HeLa cells using BCLAF1 antibody.



Western Blot analysis of various cells using BTF Polyclonal Antibody

Product Name: BTF Rabbit Polyclonal Antibody
Catalog #: APRab07685



Western Blot analysis of A549 cells using BTF Polyclonal Antibody

Note

For research use only.