

Product Name: Cdc50B Rabbit Polyclonal Antibody
Catalog #: APRab08525



Summary

Production Name	Cdc50B Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	TMEM30B
Alternative Names	TMEM30B; CDC50B; Cell cycle control protein 50B; Transmembrane protein 30B
Gene ID	161291.0
SwissProt ID	Q3MIR4.The antiserum was produced against synthesized peptide derived from human TMEM30B. AA range:91-140

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:40000. Not yet tested in other applications.
Molecular Weight	38kD

Background


Product Name: Cdc50B Rabbit Polyclonal Antibody
Catalog #: APRab08525




similarity:Belongs to the CDC50/LEM3 family.,similarity:Belongs to the CDC50/LEM3 family.,

Research Area


Image Data



Western blot analysis of lysates from K562 cells, using TMEM30B Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using TMEM30B antibody.



Western Blot analysis of K562 cells using Cdc50B Polyclonal Antibody

Product Name: Cdc50B Rabbit Polyclonal Antibody
Catalog #: APRab08525



Note

For research use only.