

Product Name: DAP-1 Rabbit Polyclonal Antibody
Catalog #: APRab09783



Summary

Production Name	DAP-1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	DAP
Alternative Names	DAP; DAP1; Death-associated protein 1; DAP-1
Gene ID	1611.0
SwissProt ID	P51397.The antiserum was produced against synthesized peptide derived from human DAP. AA range:11-60

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:40000
Molecular Weight	

Background

This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death


Product Name: DAP-1 Rabbit Polyclonal Antibody
Catalog #: APRab09783



that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014],function:Involved in mediating interferon-gamma-induced cell death.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using DAP Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.