

Product Name: DLEU7 Rabbit Polyclonal Antibody
Catalog #: APRab10014



Summary

Production Name	DLEU7 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	DLEU7 LEU7
Alternative Names	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
Gene ID	220107.0
SwissProt ID	Q6UYE1.The antiserum was produced against synthesized peptide derived from the Internal region of human DLEU7. AA range:51-100

Application

Dilution Ratio	IHC 1:50-200 ELISA 1:10000-20000
Molecular Weight	

Background

alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and


Product Name: DLEU7 Rabbit Polyclonal Antibody
Catalog #: APRab10014



kidney.,tissue specificity:Ubiquitous.,alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue specificity:Ubiquitous.,

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Note

For research use only.