

Product Name: DMRTB Rabbit Polyclonal Antibody
Catalog #: APRab10044



Summary

Production Name	DMRTB Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	DMRTB1
Alternative Names	
Gene ID	63948.0
SwissProt ID	Q96MA1.Synthesized peptide derived from human protein . at AA range: 270-350

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	37kD

Background

similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue

specificity:Testis.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue


Product Name: DMRTB Rabbit Polyclonal Antibody
Catalog #: APRab10044



specificity:Testis,.

Research Area

Image Data



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

Note

For research use only.